

See the difference **reditest**® can make.

Equip your drug test program with reliable and rapid on-site screening devices.

To place an order, call: **1-877-444-0049**. *Can't find what your looking for? Just ask us!*

With Reditest, we offer you more than the highest quality, most cost-effective devices in the industry.

We offer complete drug testing confidence.

Reditest® RediCup® URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL(S)	REDICUP® CONFIGURATIONS — 25 per box
01 102 0026	4	M-AMP/COC/THC/OPI (300)
01 102 0027	5	M-AMP/COC/THC/OPI (300)/BZO
01 102 0028	5	M-AMP/COC/THC/OPI (300)/PCP
01 102 0121	5	M-AMP/COC/THC/OPI (300)/AMP
01 102 0135	6	M-AMP/COC/THC/AMP/OPI2 (2,000)/BZO
01 102 0029	6	M-AMP/COC/THC/OPI (300)/PCP/BZO
01 102 0058	10	AMP/COC/THC/M-AMP/OPI2 (2,000)/PCP/BZO/TCA/BAR/MTD
01 102 0059	10	M-AMP/COC/THC/BZO/TCA/BAR/PCP/AMP/OPI (300)/MTD
01 102 0137	10	M-AMP/COC/THC/BZO/MTD/BAR/MDMA/OPI (300)/PCP/OXY

iCup® URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL(S)	ICUP A.D CONFIGURATIONS (with Specimen Validity Test) — 25 per box
01 102 2032	4	COC/THC/OPI (2,000)/M-AMP—CLIA Waived
01 102 2033	4	COC (150)/THC/AMP/M-AMP (500)
01 102 2034	5	COC/THC/OPI (2,000)/AMP/M-AMP—CLIA Waived
01 102 2021	5	COC/THC/MOP (300)/AMP/M-AMP
01 102 2022	6	COC/THC/OPI (2,000)/AMP/M-AMP/BZO
01 102 2037	6	COC/THC/OPI (2,000)/AMP (300)/OXY/MDMA
01 102 2038	8	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR
01 102 2069	8	COC/THC/MOP (300)/AMP/M-AMP/BZO/OXY/PCP
01 102 2039	9	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD
01 102 2129	10	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD/TCA
01 102 2074	10	COC/THC/OPI (2,000)/AMP/M-AMP/BZO/BAR/OXY/MTD/PPX
01 102 2027	12	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD/TCA/OXY/PPX
01 102 2028	13	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD/TCA/OXY/PPX/BUP
PART NUMBER	PANEL	ICUP ALL INCLUSIVE CUP CONFIGURATION
01 102 2020	10	COC/THC/OPI (2,000)/AMP/M-AMP/BZO/BAR/OXY/MDMA/PPX

E-Z Split Key Cup II URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER **PANEL** **EZ SPLIT KEY CUP II CONFIGURATION (WITH SPECIMEN VALIDITY TEST)** — 25 per box

01 102 1974 5 COC/THC/OPI (2,000)/AMP/PCP—CLIA Waived

PART NUMBER **PANEL(S)** **EZ SPLIT KEY CUP II CONFIGURATIONS** — 25 per box

01 102 1984 6 COC/THC/OPI (2,000)/AMP/M-AMP/BZO

01 102 1985 10 COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD/MDMA

01 102 2096 12 COC/THC/MOP (300)/AMP/M-AMP/BZO/BAR/MTD/MDMA/OXY/BUP/PPX

QuikScreen™ Cup URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER **PANEL** **QUICKSCREEN® CUP CONFIGURATIONS** — 25 per box

01 357 6001 11 COC/THC/OPI/AMP/M-AMP/PCP/BZO/BAR/OXY/MTD/TCA

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART # **PANEL(S)** **PANEL-DIP CONFIGURATIONS** — 25 per box

01 102 0018 1 Amphetamines (AMP)

01 102 0019 1 Barbiturates (BAR)

01 102 0173 1 Buprenorphine (BUP)

01 102 0022 1 Benzodiazepines (BZO)

01 102 0189 1 Cocaine (COC 150)

01 102 0001 1 Cocaine (COC 300)

01 102 2143 1 EDDP (Methadone Metabolite)

01 102 0190 1 Methamphetamines (M-AMP 500)

01 102 0002 1 Methamphetamines (M-AMP 1,000)

01 102 0020 1 Methadone (MTD)

01 102 0036 1 MDMA (Ecstasy)

01 102 0003 1 Opiates (OPI 300)

01 102 0037 1 Oxycodone (OXY)

01 102 0021 1 Phencyclidine (PCP)

01 102 0023 1 Tricyclic antidepressants (TCA)

01 102 0004 1 Marijuana (THC)

01 102 0007 2 COC (300)/M-AMP

01 102 0005 2 COC (300)/OPI (300)

01 102 0006 2 COC (300)/THC

continued on next page.

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART #	PANEL(S)	PANEL-DIP CONFIGURATIONS — 25 per box
01 102 0008	2	M-AMP/THC
01 102 0030	2	M-AMP/OPI (300)
01 102 0192	2	M-AMP (500)/THC (50)
01 102 0191	2	COC (150)/THC (50)
01 102 0009	3	COC (300)/M-AMP/THC
01 102 0010	3	COC (300)/OPI (300)/THC
01 102 0011	3	M-AMP/OPI (300)/THC
01 102 0014	3	COC (300)/M-AMP/OPI (300)
01 102 0193	3	COC (150)/M-AMP (500)/THC (50)
01 102 0194	3	COC (150)/OPI (300)/THC (50)
01 102 0012	4	COC (300)/M-AMP/OPI (300)/THC
01 102 0032	4	AMP (1,000)/COC (300)/OPI (300)/THC
01 102 0199	4	AMP (1,000)/COC (150)/OPI (300)/THC (50)
01 102 0195	4	COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0015	5	BZO/COC (300)/M-AMP/OPI (300)/THC
01 102 0013	5	COC (300)/M-AMP/OPI (300)/PCP/THC
01 102 0033	5	AMP/COC (300)/OPI (300)/PCP/THC
01 102 0034	5	AMP (1,000)/COC (300)/M-AMP/OPI (300)/THC
01 102 0047	5	AMP (1,000)/COC (300)/OPI (2,000)/PCP/THC
01 102 0201	5	AMP (1,000)/COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0197	5	COC (150)/BZO (300)/M-AMP (500)/OPI (300)/THC (50)
01 102 0196	5	COC (150)/M-AMP (500)/OPI (300)/PCP(25)/THC (50)
01 102 0200	5	AMP (1,000)/COC (150)/OPI (300)/PCP(25)/THC (50)
01 102 0016	6	BZO/COC (300)/M-AMP/OPI (300)/PCP/THC
01 102 0017	6	BZO/COC (300)/M-AMP/MTD/OPI (300)/THC
01 102 0024	6	BAR/BZO/COC (300)/M-AMP/OPI (300)/THC
01 102 0119	6	BZO/COC (300)/M-AMP/OPI (300)/OXY/THC
01 102 0174	6	AMP (300)/COC (150)/M-AMP (500)/MDMA/OPI (300)/THC
01 102 0175	6	BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/THC
01 102 0198	6	BZO (300)/COC (150)/M-AMP (500)/MTD(300)/OPI (300)/THC (50)
01 102 0203	6	AMP (1,000)/BZO (300)/COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0202	6	BZO (300)/COC (150)/M-AMP (500)/OPI (300)/OXY (100)/THC (50)
01 102 0035	7	AMP (1,000)/BZO/COC (300)/OPI (300)/PCP/TCA/THC
01 102 0176	7	BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/THC
01 102 0177	7	AMP/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/THC
01 102 0178	7	AMP (1,000)/COC (150)/M-AMP (500)/MDMA/OPI (300)/PCP/THC

continued on next page.

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART #	PANEL(S)	PANEL-DIP CONFIGURATIONS — 25 per box
01 102 0169	8	AMP/BZO/COC/M-AMP/MDMA/OPI (300)/OXY/THC
01 102 0179	8	AMP/BZO/COC/M-AMP/OPI (300)/OXY/PCP/THC
01 102 0180	9	AMP/BUP/BZO/COC/M-AMP/OPI (300)/OXY/PCP/THC
01 102 0181	9	AMP (300)/BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/PCP/THC
01 102 0025	10	AMP (1,000)/BAR/BZO/COC (300)/M-AMP/MTD/OPI (300)/PCP/TCA/THC
01 102 0138	10	BAR/BZO/COC (300)/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0182	10	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/THC
01 102 0183	10	BAR/BZO/COC (150)/M-AMP (500)/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0184	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/PCP/THC
01 102 0185	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (2,000)/OXY/PCP/THC
01 102 0186	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/PPX/THC
01 102 0187	11	AMP (300)/BAR/BZO/COC (150)/M-AMP (500)/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0141	12	AMP (1,000)/BAR/BZO/COC (300)/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/PPX/THC
01 102 0188	12	AMP/BAR/BUP/BZO/COC/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/THC

Reditest® 6 Cassette CLIA WAIVED URINE SUBSTANCE ABUSE SCREENING DEVICE

PART NUMBER	PANEL	REDITEST® 6 CASSETTE CONFIGURATION — 25 per box
01 102 0166	6	AMP/COC/M-AMP/OPI2 (2,000)/PCP/THC—CLIA Waived

Reditest® On-Site Oral ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICE

PART NUMBER	PANEL	REDITEST® ON-SITE ORAL CONFIGURATION — 25 per box
01 102 0127	6	AMP/COC/M-AMP/OPI/PCP/THC

iScreen® OFD ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL(S)	ISCREEN® OFD CONFIGURATIONS — 25 per box
01 102 2024	5	COC/THC/OPI/M-AMP/AMP
01 102 2025	6	COC/THC/OPI/M-AMP/AMP/PCP

OrAlert™ ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL(S)	ORALERT CONFIGURATIONS — 25 per box
01 102 1960	6	AMP/COC/M-AMP/OPI/PCP/THC
01 102 2083	6	AMP/COC/M-AMP/OPI/BZO/THC

Reditest® Alcohol-Saliva Test ORAL FLUID ALCOHOL SCREENING DEVICE

PART NUMBER	DEVICE	REDITEST® ALCOHOL-SALIVA TEST STRIP — 25 per box
01 362 0001	TEST STRIP	Detects 0.02%/ 0.08%/ AND 0.30% SAC

Reditest® Alcohol-Breath Test BREATH ALCOHOL SCREENING DEVICE

PART NUMBER	DEVICE	REDITEST® ALCOHOL-BREATH TEST DEVICE — 25 per box
01 215 0004	TEST STRIP	Detects 0.02%/ 0.08%/ AND 0.30% SAC

ALCO-SCREEN® ORAL FLUID ALCOHOL SCREENING DEVICES

PART NUMBER	DEVICE	ALCO-SCREEN® // ALCO-SCREEN® 02 — 25 per box
01 094 0055	TEST STRIP	Detects 0.02%/ 0.08%/ AND 0.30% SAC
01 094 0056	TEST STRIP	Detects 0.02%—DOT APPROVED

Reditest® Smoke Cassette URINE COTININE SCREENING DEVICES

PART NUMBER	DEVICE	REDITEST® SMOKE CASSETTE — 25 per box
01 102 0140	CASSETTE	Cotinine (200 ng/ml)—DETECTS NICOTINE METABOLITE

Adulteration Strips SPECIMEN VALIDITY TEST (SVT)

PART NUMBER	DEVICE	ADULTERATION STRIPS — 25 strips per canister
01 102 1910	TEST STRIP	1-step validity test. Screens 6 parameters for specimen validity, dilution & tampering. Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, and Oxidants/Pyridinium Chlorochromate

screening device cut-off levels

The Reditest devices detect a wide range of drugs and drug metabolites in urine or oral fluids at the following cut-off concentrations:

URINE SCREENING DEVICE CUT-OFF LEVELS

DRUG	CALIBRATOR	CUT-OFF LEVEL
Amphetamine (AMP 1,000)	d-Amphetamine	1,000 ng/mL
Amphetamine (AMP 300)	d-Amphetamine	300 ng/mL ¹
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC 300)	Benzoylcegonine	300 ng/mL
Cocaine (COC 150)	Benzoylcegonine	150 ng/mL ²
EDDP (Methadone Metabolite)	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300 ng/mL
Marijuana (THC)	11-nor- Δ 9-THC-9 COOH	50 ng/mL ²
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (M-AMP 1,000)	d-Methamphetamine	1,000 ng/mL
Methamphetamine (M-AMP 500)	d-Methamphetamine	500 ng/mL ²
Methylenedioxymethamphetamine (MDMA)	D,L Methylenedioxymethamphetamine	500 ng/mL ²
Opiates (OPI2 2,000 / OPI 2,000)	Morphine	2,000 ng/mL ²
Opiates (OPI 300 / MOP 300)	Morphine	300 ng/mL ¹
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL ¹
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL

SALIVA SCREENING DEVICE CUT-OFF LEVELS

DRUG	CALIBRATOR	CUT-OFF LEVEL
Amphetamine (AMP)	d-Amphetamine	50 ng/mL
Methamphetamine (M-AMP)	d-Methamphetamine	50 ng/mL
Benzodiazepines (BZO)	Oxazepam	10 ng/mL
Cocaine (COC)	Benzoylcegonine	20 ng/mL
Opiates (OPI)	Morphine	40 ng/mL
Marijuana (THC)	11-nor- Δ 9-THC-9 COOH	12 ng/mL ³
Marijuana (THC)	Δ 9-THC	100 ng/mL ⁴
Phencyclidine (PCP)	Phencyclidine	10 ng/mL

1. Lower than SAMHSA recommended guidelines

2. Meets SAMHSA cut-off levels

3. Reditest On-Site Oral and iScreen OFD THC detection

4. OrAlert THC detection. Includes improved sensitivity of parent THC.

5. RTL's product training certification ensures that collection, device testing and result interpretation are performed in a manner consistent with manufacturer recommendation

RediCup[®], Panel-Dip, 6 Cassette devices, EZ-Split Key Cup, iCup[®] and Alcohol-Breath Test are for professional in vitro diagnostic use only. On-Site Oral, Alcohol-Saliva Strips, iScreen[™] OFD, OrAlert[™] and ALCO-SCREEN are for forensic use only. Reditest[®] Smoke Cassette is for determination of smoking status only. ALCO-SCREEN O2 is for DOT use. The assays provide only a preliminary test result. Any positive result obtained with these drug screens are presumptive and should be confirmed by an alternate method such as GC/MS.

11 000 0135 REV19