

# A reliable and cost effective on-site cotinine screening solution.

ARE THEY SMOKING? NOW YOU'LL KNOW.

Now you can easily and quickly determine smoking status with more certainty and less cost. The Reditest® Smoke Cassette is a rapid, one step test for the qualitative detection of Cotinine (Nicotine metabolite) in human urine.

The Reditest Smoke Cassette is ideal for the workplace, corrections and schools. For the scientific confirmation of preliminary positive results, Redwood Toxicology Laboratory offers quick turnaround.

*For more information, call today!*

**877-444-0049**


[www.redwoodtoxicology.com](http://www.redwoodtoxicology.com)

## Features & benefits

- » Convenient—simple to use and easy to read
- » Quick Analysis—on-site results in 5 minutes
- » Smoking Status—detects Cotinine, a Nicotine metabolite, in urine at 200 ng/mL

## How it works.

PRODUCT TRAINING AND CERTIFICATION AVAILABLE ONLINE.



PROCEDURE NOTE: Refer to product insert for complete instructions, limitations, and warnings.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS), gas chromatography/tandem mass spectrometry (GC/MS/MS) and liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Clinical consideration and professional judgement should be applied to any drug of abuse test result, particularly when preliminary positive results are indicated.

**FOR DETERMINATION OF SMOKING STATUS ONLY.**