

See the difference **reditest**® can make.

Equip your drug test program with reliable and rapid on-site screening devices.

To place an order, call: **877.444.0049**. Can't find what you're looking for? Just ask us!

With Reditest, we offer you more than the highest quality, most cost-effective devices in the industry.

We offer complete drug testing confidence.

Alere iCup® A.D. URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL(S)	ICUP A.D. CONFIGURATIONS (with Specimen Validity Test) — 25 per box
01 102 2033	4	AMP1000/COC150/mAMP500/THC50 + S.V.T. OX,CR,PH
01 102 2032	4	COC300/mAMP1000/OPI2000/THC50 + S.V.T. OX,SG,PH
01 102 2034	5	AMP1000/COC300/mAMP1000/OPI2000/THC50 + S.V.T. OX,SG,PH
01 102 2035	5	AMP1000/COC300/OPI2000/PCP25/THC50 + S.V.T. OX,SG,PH
01 102 2021	5	AMP1000/COC300/mAMP1000/MOP300/THC50 + S.V.T. OX,SG,PH
01 102 2036	5	COC300/mAMP1000/OPI2000/PCP25/THC50 + S.V.T. OX,SG,PH
01 102 2022	6	AMP1000/BZO300/COC300/mAMP1000/OPI2000/THC50 + S.V.T. OX,SG,PH
01 102 2037	6	AMP300/COC300/MDMA500/OPI2000/OXY100/THC50 + S.V.T. OX,SG,PH
01 102 2023	6	AMP1000/COC300/mAMP1000/OPI2000/PCP25/THC50 + S.V.T. OX,SG,PH
01 102 2038	8	AMP1000/BAR300/BZO300/COC300/mAMP1000/OPI2000/PCP25/THC50 + S.V.T.OX,SG,PH
01 102 2069	8	AMP1000/BZO300/COC300/mAMP1000/MOP300/OXY100/PCP25/THC50 + S.V.T. OX,CR,PH
01 102 2291	8	MDMA500/MOP300/OXY100/BZO300/BAR300/mAmp1000/COC300/THC50 OX,CR,PH
01 102 2039	9	AMP1000/BAR300/BZO300/COC300/mAMP1000/MTD300/OPI2000/PCP25/THC50 + S.V.T. OX,SG,PH
01 102 2074	10	AMP1000/BAR300/BZO300/COC300/mAMP1000/MTD300/OPI2000/OXY100/PPX300/THC50 + S.V.T. OX,CR,PH
01 102 2129	10	AMP1000/BAR300/BZO300/COC300/mAMP1000/MTD300/OPI2000/PCP25/TCA1000/THC50 + S.V.T. OX,SG,PH,NI,GL,CR
01 102 2027	12	AMP1000/BAR300/BZO300/COC300/mAMP1000/MTD300/OPI2000/OXY100/PCP25/PPX300/TCA1000/THC50 + S.V.T. OX,CR,PH

Alere iCup® URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL	ICUP ALL INCLUSIVE CUP CONFIGURATION — 25 per box
01 102 2064	3	COC300/OPI2000 /THC50
01 102 2061	7	AMP1000/COC300/mAMP1000/MDMA500/OPI2000/PCP25/THC50
01 102 2295	9	MDMA500/MOP300/PCP25/BZO300/BAR300/COC150/MAMP500/THC50/OXY100
01 102 2292	10	MAMP1000/COC300/THC50/MDMA500/MOP300/PCP25/BZO300/BAR300/MTD300/OXY100
01 102 2293	10	MAMP1000/COC300/THC50/MDMA500/MOP300/BUP10/BZO300/BAR300/MTD300/OXY100
01 102 2020	10	AMP1000/BAR300/BZO300/COC300/mAMP1000/MDMA500/OPI2000/OXY100/PPX300/THC50
01 102 2055	10	AMP1000/BAR300/BZO300/COC300/mAMP1000/MTD300/OPI2000/PCP25/TCA1000/THC50
01 102 2294	11	MAMP1000/COC300/THC50/MDMA500/MOP300/PCP25/BZO300/BAR300/MTD300/OXY100/ BUP10
01 102 2028	13	AMP1000/BAR300/BUP10/BZO300/COC300/mAMP1000/MTD300/OPI2000/OXY100/PCP25/PPX300/TCA1000/THC50

The Alere™ Instrumented Test System with DxLINK™

PART NUMBER	PRODUCT	THE ALERE™ INSTRUMENTED TEST SYSTEM WITH DXLINK™ — FOR FORENSIC USE ONLY
01 506 0001	Scanner	Technology combines a novel software system with an image capture device allowing for automated test interpretation of onsite forensic iCup drug test devices. For compatible devices, see iCups listed below.
PART NUMBER	PANEL(S)	ICUP ALL INCLUSIVE CUP CONFIGURATION — 25 per box
01 102 2283	8	AMP1000/BZO300/COC300/mAMP1000/MOP300/OXY100/PCP25/THC50 + S.V.T. OX,CR,PH— For Forensic Use Only
01 102 2281	10	BAR300/BUP10/BAR300/COC300/MAMP1000/MDMA500/MOP300/MTD300/OXY100/THC50— For Forensic Use Only
01 102 2284	13	AMP1000/BAR300/BUP10/BZO300/COC300/mAMP1000/MTD300/OPI2000/OXY100/PCP25/PPX300/TCA1000/THC50— For Forensic Use Only

E-Z Split Key Cup II URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL	EZ SPLIT KEY CUP II CONFIGURATION (WITH SPECIMEN VALIDITY TEST) — 25 per box
01 102 1974	5	COC/THC/OPI (2,000)/AMP/PCP—CLIA Waived
PART NUMBER	PANEL(S)	EZ SPLIT KEY CUP II CONFIGURATIONS — 25 per box
01 102 1984	6	COC/THC/OPI (2,000)/AMP/M-AMP/BZO
01 102 1985	10	COC/THC/OPI (2,000)/AMP/M-AMP/PCP/BZO/BAR/MTD/MDMA
01 102 2096	12	COC/THC/MOP (300)/AMP/M-AMP/BZO/BAR/MTD/MDMA/OXY/BUP/PPX

QuikScreen™ Cup URINE SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER	PANEL	QUICKSCREEN® CUP CONFIGURATIONS — 25 per box
01 357 6001	11	COC/THC/OPI/AMP/M-AMP/PCP/BZO/BAR/OXY/MTD/TCA

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART #	PANEL(S)	PANEL-DIP CONFIGURATIONS — 25 per box
01 102 0018	1	Amphetamines (AMP)
01 102 0019	1	Barbiturates (BAR)
01 102 0173	1	Buprenorphine (BUP)
01 102 0022	1	Benzodiazepines (BZO)
01 102 0189	1	Cocaine (COC 150)
01 102 0001	1	Cocaine (COC 300)
01 102 0190	1	Methamphetamines (M-AMP 500)
01 102 0002	1	Methamphetamines (M-AMP 1,000)
01 102 0020	1	Methadone (MTD)
01 102 0036	1	MDMA (Ecstasy)

continued on next page.

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART #PANEL(S)PANEL-DIP CONFIGURATIONS — 25 per box

01 102 0003	1	Opiates (OPI 300)
01 102 0037	1	Oxycodone (OXY)
01 102 0021	1	Phencyclidine (PCP)
01 191 6335	1	Synthetic Cannabinoids (SynCANN)— For Forensic Use Only
01 102 0023	1	Tricyclic antidepressants (TCA)
01 102 0004	1	Marijuana (THC)
01 102 0007	2	COC (300)/M-AMP
01 102 0005	2	COC (300)/OPI (300)
01 102 0006	2	COC (300)/THC
01 102 0008	2	M-AMP/THC
01 102 0030	2	M-AMP/OPI (300)
01 102 0192	2	M-AMP (500)/THC (50)
01 102 0191	2	COC (150)/THC (50)
01 102 0009	3	COC (300)/M-AMP/THC
01 102 0010	3	COC (300)/OPI (300)/THC
01 102 0011	3	M-AMP/OPI (300)/THC
01 102 0014	3	COC (300)/M-AMP/OPI (300)
01 102 0193	3	COC (150)/M-AMP (500)/THC (50)
01 102 0194	3	COC (150)/OPI (300)/THC (50)
01 102 0012	4	COC (300)/M-AMP/OPI (300)/THC
01 102 0032	4	AMP (1,000)/COC (300)/OPI (300)/THC
01 102 0199	4	AMP (1,000)/COC (150)/OPI (300)/THC (50)
01 102 0195	4	COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0015	5	BZO/COC (300)/M-AMP/OPI (300)/THC
01 102 0013	5	COC (300)/M-AMP/OPI (300)/PCP/THC
01 102 0033	5	AMP/COC (300)/OPI (300)/PCP/THC
01 102 0034	5	AMP (1,000)/COC (300)/M-AMP/OPI (300)/THC
01 102 0047	5	AMP (1,000)/COC (300)/OPI (2,000)/PCP/THC
01 102 0201	5	AMP (1,000)/COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0197	5	COC (150)/BZO (300)/M-AMP (500)/OPI (300)/THC (50)
01 102 0196	5	COC (150)/M-AMP (500)/OPI (300)/PCP(25)/THC (50)
01 102 0200	5	AMP (1,000)/COC (150)/OPI (300)/PCP(25)/THC (50)
01 102 0016	6	BZO/COC (300)/M-AMP/OPI (300)/PCP/THC
01 102 0017	6	BZO/COC (300)/M-AMP/MTD/OPI (300)/THC
01 102 0024	6	BAR/BZO/COC (300)/M-AMP/OPI (300)/THC
01 102 0119	6	BZO/COC (300)/M-AMP/OPI (300)/OXY/THC

continued on next page.

Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PART #PANEL(S)PANEL-DIP CONFIGURATIONS — 25 per box

01 102 0174	6	AMP (300)/COC (150)/M-AMP (500)/MDMA/OPI (300)/THC
01 102 0175	6	BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/THC
01 102 0203	6	AMP (1,000)/BZO (300)/COC (150)/M-AMP (500)/OPI (300)/THC (50)
01 102 0202	6	BZO (300)/COC (150)/M-AMP (500)/OPI (300)/OXY (100)/THC (50)
01 102 0035	7	AMP (1,000)/BZO/COC (300)/OPI (300)/PCP/TCA/THC
01 102 0176	7	BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/THC
01 102 0177	7	AMP/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/THC
01 102 0169	8	AMP/BZO/COC/M-AMP/MDMA/OPI (300)/OXY/THC
01 102 0179	8	AMP/BZO/COC/M-AMP/OPI (300)/OXY/PCP/THC
01 102 0180	9	AMP/BUP/BZO/COC/M-AMP/OPI (300)/OXY/PCP/THC
01 102 0181	9	AMP (300)/BZO/COC (150)/M-AMP (500)/MDMA/OPI (300)/OXY/PCP/THC
01 102 0025	10	AMP (1,000)/BAR/BZO/COC (300)/M-AMP/MTD/OPI (300)/PCP/TCA/THC
01 102 0138	10	BAR/BZO/COC (300)/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0182	10	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/THC
01 102 0183	10	BAR/BZO/COC (150)/M-AMP (500)/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0184	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/PCP/THC
01 102 0185	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (2,000)/OXY/PCP/THC
01 102 0186	11	AMP/BAR/BUP/BZO/COC/M-AMP/MTD/OPI (300)/OXY/PPX/THC
01 102 0187	11	AMP (300)/BAR/BZO/COC (150)/M-AMP (500)/MDMA/MTD/OPI (300)/OXY/PCP/THC
01 102 0141	12	AMP (1,000)/BAR/BZO/COC (300)/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/PPX/THC
01 102 0188	12	AMP/BAR/BUP/BZO/COC/M-AMP/MDMA/MTD/OPI (300)/OXY/PCP/THC

Reditest® 6 Cassette CLIA WAIVED URINE SUBSTANCE ABUSE SCREENING DEVICE

PART NUMBER PANEL REDITEST® 6 CASSETTE CONFIGURATION — 25 per box

01 102 0166	6	AMP/COC/M-AMP/OPI2 (2,000)/PCP/THC—CLIA Waived
-------------	---	--

Alere iScreen® OFD ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER PANEL(S) ISCREEN® OFD CONFIGURATIONS — FOR FORENSIC USE ONLY — 25 per box

01 102 2024	5	COC/THC/OPI/M-AMP/AMP
01 102 2025	6	COC/THC/OPI/M-AMP/AMP/PCP

OrAlert® ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PART NUMBER PANEL(S) ORALERT CONFIGURATIONS — FOR FORENSIC USE ONLY — 25 per box

01 102 1960	6	AMP/COC/M-AMP/OPI/PCP/THC
01 102 2083	6	AMP/COC/M-AMP/OPI/BZO/THC

AlcoMate® Premium 7000 BREATH ALCOHOL SCREENING DEVICE

PART NUMBER	PRODUCT	DESCRIPTION
01 581 7000	AlcoMate® Premium Kit DOT Approved	Kit includes: Breathalyzer, User's Manual, 2 x AA Batteries, 55 x mouth pieces, 1 pre-installed Sensor Module (PRISM™), 1 additional Sensor Module, hard plastic carrying case, car adapter, handstrap
01 581 7001	AlcoMate® Sensor Module	1 Pre-calibrated replacement Sensor Module for AlcoMate Premium device
01 581 7050	AlcoMate® Mouthpieces	50-Pack of standard mouthpieces for AlcoMate Premium device

Mission® Breath Alcohol Detector BREATH ALCOHOL SCREENING DEVICE

PART NUMBER	DEVICE	MISSION® BREATH ALCOHOL DETECTOR — 510k CLEARED — 20 per box
01 532 0020	BREATH TUBE	Detects 0.02% SAC

Reditest® Alcohol-Saliva Test ORAL FLUID ALCOHOL SCREENING DEVICE

PART NUMBER	DEVICE	REDITEST® ALCOHOL-SALIVA TEST STRIP — FOR FORENSIC USE ONLY — 25 per box
01 362 0001	TEST STRIP	Detects 0.02%/ 0.08%/ AND 0.30% SAC

ALCO-SCREEN® ORAL FLUID ALCOHOL SCREENING DEVICES

PART NUMBER	DEVICE	ALCO-SCREEN® // ALCO-SCREEN® 02 — 510k CLEARED, CLIA WAIVED — 24 per box.
01 094 0055	TEST STRIP	Detects SAC at 0.02%, 0.08% and 0.03%—CLIA Waived
01 094 0056	TEST STRIP	Detects 0.02%—DOT approved, CLIA Waived

Cotinine Screening URINE & ORAL FLUID COTININE SCREENING DEVICES

PART NUMBER	DEVICE	REDITEST® SMOKE CASSETTE — 25 per box
01 102 0140	CASSETTE	Cotinine (200 ng/mL)—Detects nicotine metabolite in urine

Adulteration Strips SPECIMEN VALIDITY TEST (SVT)

PART NUMBER	DEVICE	ADULTERATION STRIPS — 25 strips per canister
01 102 1910	TEST STRIP	1-step validity test. Screens 6 parameters for specimen validity, dilution & tampering. Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, and Oxidants/Pyridinium Chlorochromate

Screening device cutoff levels

The Reditest devices detect a wide range of drugs and drug metabolites in urine or oral fluids at the following cutoff concentrations:

URINE DRUG SCREENING DEVICE CUTOFF LEVELS

DRUG	CALIBRATOR	CUTOFF LEVEL
Amphetamine (AMP 1,000)	d-Amphetamine	1,000 ng/mL
Amphetamine (AMP 300)	d-Amphetamine	300 ng/mL ¹
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC 300)	Benzoyllecgonine	300 ng/mL
Cocaine (COC 150)	Benzoyllecgonine	150 ng/mL ²
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	50 ng/mL ²
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (M-AMP 1,000)	d-Methamphetamine	1,000 ng/mL
Methamphetamine (M-AMP 500)	d-Methamphetamine	500 ng/mL ²
Methylenedioxymethamphetamine (MDMA)	D,L Methylenedioxymethamphetamine	500 ng/mL ²
Nicotine (COT)	Cotinine	200 ng/mL
Opiates (OPI2 2,000 / OPI 2,000)	Morphine	2,000 ng/mL ²
Opiates (OPI 300 / MOP 300)	Morphine	300 ng/mL ¹
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	25 ng/mL ²
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1,000 ng/mL

ORAL FLUID DRUG SCREENING DEVICE CUTOFF LEVELS

DRUG	CALIBRATOR	CUTOFF LEVEL
Alcohol (ALC)	Ethanol	0.02%/0.08%/0.30% SAC
Amphetamine (AMP)	d-Amphetamine	50 ng/mL
Methamphetamine (M-AMP)	d-Methamphetamine	50 ng/mL
Benzodiazepines (BZO)	Oxazepam	10 ng/mL
Cocaine (COC)	Benzoyllecgonine	20 ng/mL
Opiates (OPI)	Morphine	40 ng/mL
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	12 ng/mL ³
Marijuana (THC)	Δ^9 -THC	100 ng/mL ⁴
Nicotine (COT)	Cotinine	30 ng/mL
Phencyclidine (PCP)	Phencyclidine	10 ng/mL

1. Lower than SAMHSA recommended guidelines

2. Meets SAMHSA cutoff levels

3. iScreen OFD THC detection

4. OrAlert THC detection. Includes improved sensitivity of parent THC.

Reditest® Panel-Dip, Reditest® 6 Cassette devices, EZ-Split Key Cup, Alere iCup® and Mission® Breath Alcohol Detector are for professional in vitro diagnostic use only. Alcohol-Saliva Strips, Alere iScreen™ OFD, OrAlert™ and ALCO-SCREEN are for forensic use only. Reditest® Smoke Cassette is for determination of smoking status only. ALCO-SCREEN O2 and AlcoMate® Premium 7000 Breath Alcohol Test are DOT approved. The assays provide only a preliminary test result. Any positive result obtained with these drug screens are presumptive and should be confirmed by an alternate method such as GC/MS or LC/MS/MS.